



**LEATHER
UPHOLSTERY!**



**LIFT
position**



**RECLINE
position**



BROWN
(LC0504)



TAUPE
(LC0505)



BLACK
(LC0506)

MONTREAL WALL HUGGER LIFT CHAIR

LC0504 - BROWN | LC0505 - TAUPE | LC0506 - BLACK

FEATURES

- Lift / Rest / Recline armchair
- Allows the user to be easily raised to a standing position and offers the comfort of a backrest recline
- User friendly handset to adjust chair as required
- Single motor lift
- Added comfort of a pocket spring seat base
- Luxurious leather coverings
- Available in modern colours: Brown / Taupe / Black
- Quality actuators
- Warranty:
Frame: 3 years / Electrics: 1 year / Leather: 1 year

SPECIFICATIONS

All dimensions shown are nominal

- Overall: **845mm (W) x 820mm (D) x 1145mm (H)**
Seat: **490mm (W) x 520mm (D) x 500mm (H)**
Armrest: **160mm (H) x 500mm (L)**
Backrest: **695mm (W) x 730mm (H): 93 - 125°**
Legrest: **570mm (W) x 480mm (L): 90 - 180°**
Distance from wall: **100mm**
Product Weight: **55kg**
Maximum User Weight: **150kg**

Please contact **REDGUM Brand** if you require further information

P: 08 9248 4180

E: orders@for-de.com.au

REDGUM brand - 1 Business way, Malaga, Western Australia 6090

redgumbrand.com.au

IMPORTANT NOTE:

It is the responsibility of the user to ensure the suitability of this product for its intended use.
The user should consult their appropriate healthcare professional if unsure of the correct product usage for their unique circumstances.