

Pearl Deluxe Aluminium Wheelchair

Part Number RG18AAD / RG20AAD / RG18AADRB / RG20AADRB



Features

- Aluminium anti corrosive folding frame
- Padded adjustable height arm rests
- Adjustable height rear handles
- Adjustable comfort back rest (Velcro)
- 24" Solid tyres
- 7" Front rubber tyres
- Adjustable castor angle
- 4 stage seat adjustment
- Hard wearing nylon upholstery
- Height & angle adjustable foot plates
- Easy to operate lever swing away leg rests
- Manual brake on rear wheels
- Quick disconnect rear wheels
- Anti tip wheels fitted as standard
- Available in Pearlescent White

Dimensions

- Seat Width 18" (470mm) 20" (510mm)
- Seat Depth 440mm
- Seat Height 410mm
- Max User Capacity 18" - 120kg
- Max User Capacity 20" - 130kg
- Wheelchair Weight 18" - 17.0kg
- Wheelchair Weight 20" - 18.5kg

Rear Hand Brake Model
now available

Phone: 08 9248 4180
Fax: 08 9248 9452
E-mail: orders@redgumbrand.com.au