

Hospital Wheelchair with Oxygen Bottle Holder & IV Drip Holder



Features

- Strong steel frame
- Wheelchair frame can be secured for anti-theft
- Easy clean polyurethane arm rests
- 24" Solid tyres
- 7" Front rubber tyres
- Adjustable castor angle
- Hard wearing nylon upholstery
- Anti-theft swing away foot plates
- Easy to operate lever swing away leg rests
- Manual brake on rear wheels
- Anti tip wheels fitted as standard
- Oxygen bottle holder fitted as standard
- I.V. Pole fitted as standard
- Available in Gun Metal Grey

Dimensions

- Seat Width 18" (480mm)
- Seat Depth 410mm to 445mm
- Seat Height 490mm
- Max User Capacity - 150kg
- Wheelchair Weight - 17.0kg

Phone: 08 9248 4180
Fax: 08 9248 9452
E-mail: orders@redgumbrand.com.au