

Deluxe Folding Transporter



Features

- Lightweight aluminium construction
- Powder coated folding frame
- Padded arm rests (reversible)
- Fold down backrest
- 12" Solid tyres
- 8" Front rubber tyres
- Hard wearing nylon upholstery
- Swing away leg rests
- Attendant rear control handbrakes
- Manual brake on rear wheels
- Seat belt fitted
- Rear storage pocket on backrest
- Available in marble finish Red or Blue

Dimensions

- Seat Width 17" (420mm) RG17HBD
- Seat Width 19" (470mm) RG19HBD
- Seat Width 22" (550mm) RG22HBD
- Seat Depth 390mm
- Seat Height 460mm
- Max User Capacity 17/19" - 110kg
- Max User Capacity 22" - 125kg
- Wheelchair Weight 17/19" - 10.5kg
- Wheelchair Weight 22" - 13.0kg



Phone: 08 9248 4180
Fax: 08 9248 9452
E-mail: orders@redgumbrand.com.au