

TOP SELLING PRODUCT !



REDGUM brand strives to provide the highest quality products possible that are subject to change without notice.

IMPORTANT NOTE:

It is the responsibility of the user to ensure the suitability of this product for its intended use. The user should consult their appropriate healthcare professional if unsure of the correct product usage for their unique circumstances.

KNEE WALKER

RG30KW

FEATURES

- Easy to assemble
- Padded contoured knee rest - height adjustable
- Offset knee rest position - can be left or right
- Height adjustable handlebars
- Anti slide hand grips
- Left + Right hand dual locking brake system
- Large wheels for excellent mobility + stability
- Strong lightweight construction
- Folds down for storage and transport
- Removable wire carry basket
- Walking stick holders

SPECIFICATIONS (All dimensions shown are nominal)

Overall Dimensions (at lowest setting):
830mm (L) x 450mm (W) x 790mm (H)

Knee Platform Size:
370mm (L) x 185mm (W)

Knee Platform Height:
400 - 515mm (from floor to knee pad)

Handlebar Width:
450mm

Handlebar Height:
790 - 950mm (from floor to top of handlebar)

Wheel Diameter:
200mm / 8"

Product Weight:
12kg

Maximum User Weight:
130kg

Please contact **REDGUM** if you require any more information on this product.

P: 08 9248 4180

F: 08 9248 9452

E: orders@for-de.com.au

REDGUM brand - 1 Business way, Malaga, Western Australia 6090

redgumbrand.com.au