

## Forearm Walker

RG4310AA



### Features

- Ergonomically designed
- Strong steel construction
- Powder coated finish
- Easy to Fold frame
- Padded guttered armrests
- 6 Positions for height adjustment
- Adjustable brake lever depth
- Locking brake mechanism
- 200mm micro cellular wheels
- Seat platform 600mm from floor
- Wire plastic coated basket supplied
- Clear acrylic tray supplied
- Available in marble finish Red

### Dimensions

- Armrest Height 975mm to 1115mm
- Seat / Platform Height 600mm
- Seat / Platform Width 375mm
- Unit Width 580mm
- Unit Depth 680mm
- Unit Weight 14kg
- Folded Width 580mm
- Folded Depth 300mm
- Folded Height 1210mm
- Maximum User Capacity 120kg



Phone: 08 9248 4180  
Fax: 08 9248 9452  
E-mail: [orders@redgumbrand.com.au](mailto:orders@redgumbrand.com.au)