

REDGUM BRAND Daily Living and Lifestyle Products

“CrossOver” Seat Walker & Occasional Transit Chair Walker

RG4408



Can be used with
or without
removable
leg rests



Features

- Side folding seat walker
- Converts to a transit wheelchair
- Aluminium construction
- Top quality powder coated finish
- Side folding action
- Soft operation hand brake
- Seat height at 560mm from floor
- 115mm wide padded backrest
- Removable swing away leg rests fitted
- Walking stock holder fitted as standard
- 200mm spoked wheels
- Available in three colours
- Blue / Red / / White

Dimensions

- | | |
|--------------------------------|-----------------|
| • Seat Height | 560mm |
| • Seat Width | 450mm |
| • Seat Depth | 310mm |
| • Unit Overall Width | 590mm |
| • Unit Overall Length | 980mm |
| • Handle Height Range | 900mm to 1000mm |
| • Unit Net Weight | 9.42kg |
| • Unit Net Weight no footrests | 8.3kg |
| • Wheel Size | 200mm / 8” |
| • Maximum User Capacity | 120kg |

Phone: 08 9248 4180
Fax: 08 9248 9452
E-mail: orders@redgumbrand.com.au

REDGUM BRAND
Daily Living and
Lifestyle Products

REDGUM Brand is a 100% Australian Owned & Operated