

Aluminium Side Folding Walker

RG4411



Features

- Ergonomically designed
- Lightweight aluminium construction
- Top quality powder coated finish
- Single hand folding design
- Side folding action for easy storage
- Soft operation hand brake
- Seat height at 550mm from floor
- 115mm wide padded backrest
- Walking stock holder fitted as standard
- Removable storage bag fitted as standard
- 200mm spoked wheels
- Available in four colours
- Blue / Red / Purple /White

Dimensions

- Seat Height 550mm
- Seat Width 480mm
- Unit Overall Width 615mm
- Unit Overall Length 770mm
- Handle Height Range 870mm to 970mm
- Unit Net Weight 7.0kg
- Wheel Size 200mm / 8"
- Maximum User Capacity 136kg



Phone: 08 9248 4180

Fax: 08 9248 9452

E-mail: orders@redgumbrand.com.au