

Hematite Lightweight Transit Wheelchair

RG4600



Features

- Lightweight aluminium construction
- Super Lightweight 8.0kg Net Weight
- Powder coated folding frame
- Padded arm rests
- Fold down backrest
- 8" Solid rear tyres
- 8" Front rubber tyres
- Hard wearing nylon upholstery
- Swing away leg rests
- Non skid foot plates with heel strap
- Manual brake on rear wheels
- Seat belt fitted
- Rear storage pocket on backrest
- Available in marble finish Red



Dimensions

- Seat Width 18" (460mm)
- Seat Depth 400mm
- Seat Height 480mm
- Max User Capacity - 110kg
- Wheelchair Weight - 8.0kg

Options

- Elevating Legs Rest

Phone: 08 9248 4180
Fax: 08 9248 9452
E-mail: orders@redgumbrand.com.au