

Compact Wheelchair Carrier RGCWCC



REDGUM Brand has recently updated the wheelchair carrier to the new **Compact Wheelchair Carrier**, it has been designed to be lighter and overall more compact than the previous version. Designed to carry self propelled folding wheelchairs on the back of vehicles and the wheelchair carrier can be easily folded and removed from the rear of the vehicle when not required.

A drop in basket to allow for transit wheelchairs is now also available as an option

Part number - RGCWCCB

The REDGUM Compact Wheelchair Carrier has been manufactured from durable steel for long life and is also coated to prevent corrosion.

Please Note

The wheelchair requires a Reece type hitch to accommodate the fitting to a vehicle.

Dimensions

- Max Capacity 32kg
 - Total Weight 10.5kg
- Part Number RGCWCC

Features

- Easy fold mechanism
 - Easy remove feature
 - Steel construction
- Transit Chair basket available as an option
(P/N RGCWCCB)

Phone: 08 9248 4180
Fax : 08 9248 9452
E-mail: orders@redgumbrand.com.au