

SCOOTER XTRAS

NEW PRODUCT RELEASE !



*Bottle not included

IMPORTANT NOTE:

It is the responsibility of the user to ensure the suitability of this product for its intended use.

The user should consult their appropriate healthcare professional if unsure of the correct product usage for their unique circumstances.

MOBILITY BOTTLE HOLDER

SE-020 (Black)

FEATURES

- Adjusts to take most bottle sizes
- Easy to attach / adjust and re-position
- Swivel head for angle adjustment
- Adjustable mount can attach onto any tube / ledge up to 40mm thick
- Strong / secure grip
- Strong lightweight construction

SPECIFICATIONS (All dimensions shown are nominal)

Overall Dimensions:

Width: **65mm**

Depth: **155mm**

Height: **90mm**

Bottle Size:

Ø55 - 70mm

Adjustable Grip:

0 - 40mm (will attach to round or flat surface)

Product Weight:

160g

Maximum load:

1kg / 1 litre

Product Materials:

Contact Scooter Xtras for details

Please contact **Scooter Xtras** if you require further information on this product.

P: 08 9248 4180

F: 08 9248 9458

E: orders@shoprider.com.au

Scooter Xtras - 1 Business way, Malaga, Western Australia

scooterxtras.com.au