

Titan Heavy Duty Wheelchair

RG22HD



Features

- 200kg maximum weight capacity
- Powder coated steel folding frame
- Double cross bar seat support
- Adjustable back rest (4 x 50mm straps)
- Arm rests with anti deforming pads
- Lift up armrests
- 24" Hyper strength solid tyres
- Easy remove rear wheels
- 8" Front rubber tyres
- Attendant control rear drum brakes
- Manual brake on rear wheels
- Hard wearing upholstery inc reflect tape
- East to operate swing away leg rests
- Foot plate joiner
- Anti tip wheels fitted as standard
- Available in Black finish
- Optional Stump support to suit
- Optional Elevating leg rests to suit



Dimensions

- Seat Width 22" (590mm)
- Seat Depth 430mm
- Seat Height 450mm
- Max User Capacity - 200kg
- Wheelchair Weight - 25.0kg

Phone: 08 9248 4180
Fax: 08 9248 9452
E-mail: orders@redgumbrand.com.au