



BURGUNDY RED

COBALT BLUE

889ASN

## FEATURES

- High back contoured seat with headrest
- Gas adjustable steering tiller
- Loop handle with finger operated controls
- High / Low speed settings
- Rear wheel drive
- Electromagnetic brake
- Front + Rear (Adjustable) suspension
- Twin Headlights (LED)
- Rear lights
- Indicators / Hazard lights
- Large battery capacity
- Battery power indicator
- USB power outlet
- Sliding / Rotating / Height adjustable seat
- Lift up & adjustable width arm rests
- Removable shopping basket

## SPECIFICATIONS

Overall: 1390mm (L) x 610mm (W) x 1310mm (H)

Turning Circle: 1750mm

Speed: 10km/h

Maximum range: Up to 32km

Safe climb angle: 9°

Maximum climb angle: 9°

Motor: 1.5Hp high torque motor (4 Pole)

Batteries: 2 x 50 Amp minimum (not supplied)

Charger: 8 Amp

Tyres (Pneumatic): Ø330mm - Front + Rear

Product Weight: 104kg - with batteries

Maximum User Weight: 159kg - on flat ground

